

# SPECIFICATIONS



## PLATES/LEATHERS

SIZE	A	B	C	D	E
SMALL	61	67.5	29	38.5	66.5
MEDIUM	64.5	71.5	30.5	40	71.5
LARGE	68	75	31.5	43	75
EXTRA LARGE*	68	80	31.5	47.5	81
EXTRA LONG**	-	92.5	-	-	77

\*Extra Large Plates Set comprises Large Plates with Extra Large Leather.  
 \*\*Extra Long Leather - sold as Leather Only

## SPACERS


SIZE	17	18	19	20	21	22	23	24	25
mm	56.5	59.5	62.5	66	69	72.5	76	79	82

ALL CHART MEASUREMENTS SHOWN ARE IN MILLIMETRES

# SELECTING YOUR SPACER SIZE

## GETTING STARTED

You'll need a 10cm wide strip of paper, a pen or fine tip marker, and a ruler.  
 First take a look at your middle finger - the one which will go through the ring. You need to measure the part of this finger where the ring will sit, just behind the knuckle



## MEASURE

Now, take your strip of paper and wrap it around the section of your finger that you need to measure. Take your pen/marker and draw across the paper strip at the point where the end of the paper overlaps, measuring the circumference of your finger.



## COMPARE

Straighten out the strip and measure the length between the marks where the overlaps were. This length now can be compared with the chart above to select your correct spacer size. Generally if you are between sizes then go for the smaller size.

I strongly recommend you do this once yourself, then get someone to double check your measurements.

\*The ring is made of slightly stretchy rubber, however if you do have larger knuckles you may need to go up a size.

